


SCANDITRONIX

WELLHÖFER

A smart modular phantom

ImRT Phantom



Member of the IBA Group

The phantom solution for IMRT

FAST_{est} | **most ACCURATE** | **most RELIABLE**

- ✓ Modular design for maximum flexibility
- ✓ Universal Body, Head and Neck as well as Stereotactic verification
- ✓ Film, Absolute Dosimetry and TLD

- ✓ Multiple film exposure
- ✓ Consists of water equivalent RW3 material (1.045g/m³ density)
- ✓ Easy and fast adjustment under Linac and in CT scanners


- ✓ Developed in scientific cooperation with the University Hospital Hamburg-Eppendorf
- ✓ Proven clinical results with Royal Marsden Hospital (Sutton/United Kingdom), Institut Gustave Roussy (Villejuif/France) and Ospedale S. Anna (Como/Italy)

During implementation of IMRT treatments in clinical routine, film dosimetry is used to verify the planned TPS dose versus the measured dose distribution.

The I'mRT Phantom is a water equivalent phantom for multiple film measurements and verification of the absolute dose. The unique modular design allows Universal Body as well as Head and Neck and Stereotactic applications.

Universal Body verification


The I'mRT Phantom consists of a universal body shaped section in which up to 15 films can be exposed simultaneously to larger IMRT fields. Subsequent comparison of the film dose distribution with the planned TPS data in the OmniPro-I'mRT software provides an accurate verification of the IMRT delivery process.



Head and Neck/Stereotactic verification

The modular *ImRT* Phantom design also includes a removable cubic phantom that has been developed in cooperation with the University Hospital Hamburg-Eppendorf and with Euromechanics GmbH (Schwarzenbruck/Germany). It is specifically designed for Head and Neck and Stereotactic applications. Markers in different colors are engraved on the surface of the cube for easy adjustment of the phantom under the Linac and in a CT scanner.

With film distance plates of 1 cm thickness up to 15 films with a maximum size of 16 cm x 16 cm can be positioned in a transversal, coronal or sagittal orientation. The films are cut in the darkroom, alternatively radiochromic films can be used. The cubic phantom docks into the torso, thus restoring the original anthropomorphic shape.



Absolute Dose verification

For absolute dose verification together with DOSE-1 different types of ionization chambers can be inserted at all points in the cubic phantom module within a 1 cm grid. Several auxiliary inserts with different depths and for different chamber types are available. In combination with the film distance plates the ionization chambers can be placed at any position in the phantom.




Options

Optionally, a set of 6 localizer plates is available for the use of the cubic phantom in a CT scanner. As another option one or more TLD plates can be inserted in the cubic phantom, either as stand-alone application or in conjunction with film. Each TLD plate can be loaded with up to 196 TLD detector rods on each side with an individual spacing of 1 cm. Furthermore TLD cubes and other TLD inserts are available on request.



CT localizer plates



TLD plate

Technical specifications

ImRT Phantom

ImRT Phantom

Size of complete phantom:	33 cm (L) x 36 cm (W) x 18 cm (H)
Carriage and levelling plate:	33 cm (L) x 44 cm (W) x 1 cm (H)
Material:	RW3
Density:	1.045 g/cm ³
Weight:	22 kg

Universal Body section

Number of slabs:	15
Thickness:	1 cm
Film type:	ready-packed films/radiochromic films
Number of films:	up to 15
Film spacing:	min. 1 cm
Geometry:	transversal
Marking of films:	universal Body film marker. 3 marking holes define phantom origin, axes and film orientation

Head and Neck/Stereotactic cubic section

Dimensions (outer):	18 cm (L) x 18 cm (W) x 18 cm (H)
Film size:	16 cm (L) x 16 cm (W)
Film spacing:	min. 1 cm
Number of films:	up to 15
Geometry:	transversal, coronal or sagittal orientation
Compensation of film thickness:	distance plates included: 16 plates with 10 mm thickness. 1 plate with 1 mm thickness, 2 plates with 2 mm thickness and 1 plate with 5 mm thickness
Film preparation:	radiographic film cut in darkroom or radiochromic film
Film marking:	manually
Ionization chamber inserts:	Farmer type (FC65G/P) included
Chamber positioning:	at all points in the cube within a 1 cm grid
Lateral scattering bodies:	two lateral scattering bodies that can be mounted on all sides of the cube
Width with scattering bodies:	36 cm

Options

Ionization chamber inserts:	Compact type and inserts for detectors of other manufacturers available on request
CT localizer plates:	set of 6 localizer plates for the use of the cubic phantom in a CT scanner
TLD plate inserts:	up to 196 TLD detector rods (Ø 1 * 6 mm), 1 cm spacing/TLD cubes/other TLD inserts on request

Technical data is subject to change without prior notice.

Sweden: Scanditronix Wellhöfer AB – Stålgatan 14 – SE-754 50 Uppsala
Tel.: +46 18 18 07 00 ☐ Fax: +46 18 12 75 52

Germany: Scanditronix Wellhöfer GmbH – Bahnhofstr. 5 – DE-90592 Schwarzenbruck
Tel.: +49 9128 607 0 ☐ Fax: +49 9128 607 10

USA: Scanditronix Wellhofer North America – 3150 Stage Post Drive, Suite 110 – Bartlett, IN 38133
Tel.: +1 901 386 2242 ☐ Fax: +1 901 382 9453

China: Scanditronix Wellhöfer China – 100015-92 Dashanzi, Chaoyang district – CN-100015 Beijing
Tel.: +86 10 6437 2760 ☐ Fax: +86 10 6437 2763

www.scanditronix-wellhofer.com ☐ info@scanditronix-wellhofer.com